


| | | |
|---|---|---------|
|  DESSILLONS & DUTRILLAUX - Since 1949 | TECHNICAL CHARACTERISTIC | CT 036 |
| | Non-invasive blood pressure cuff made in "Soft Cuff®" fabric | Index 7 |

Regulatory information

| | | |
|------------|-------------------------|--|
| CE Marking | Index of classification | <div style="border: 1px solid black; padding: 5px; display: inline-block;"> PHTALATES LATEX </div> |
| CE | I | |



The range of cuffs Soft Cuff is intended to be used to measure non invasive blood pressure associated with electronic or manual blood pressure.

This flexible and robust design can be cleaned in a washing machine to clean them by immersion with the requirement to remove the pocket.

- Reliability : an Index and a range "Min-Max" can quickly validate if the product is adapted to the size of the patient.
 The indication "artery" can perfectly position the cuff, pocket located automatically centered on the passage of the artery.

- Safety : to fight against nosocomial infections, fabric cuff is antimicrobial with a specific treatment.

1 - Device identification


Blood pressure cuff 1 tube



Blood pressure cuff 2 tubes




| | | | |
|---|---|--|------------|
| Dessillons & Dutrillaux Zone Industrielle La Tuque Tél : +33 5 53 48 30 66 technique@ddmedical.fr | 47240 CASTELCULIER Fax : +33 5 53 47 24 64 www.ddmedical.fr | Date of first EC marking certificate : 2004 | Page 1 / 8 |
|---|---|--|------------|

| | | |
|---|---|---------|
|  DESSILLONS & DUTRILLAUX - Since 1949 | TECHNICAL CHARACTERISTIC | CT 036 |
| | Non-invasive blood pressure cuff made in "Soft Cuff®" fabric | Index 7 |








2 – References

According to EN 1060-1 EC, the optimal size of the inflation bag must be at least equal in width to 40% of the circumference of the arm, and a length between 80 and 100% of the circumference of the arm (see markings shown on cuffs Soft Cuff).


| Reference | Designation | Tube | Circumference in cm | | Dimensions in mm | | | | | Tube color |
|-----------|-----------------------|------|---------------------|-----|------------------|-----|-----|-----|-----|------------|
| | | | Min | Max | A | B | C | D | E | |
| BT101 | New born | 1 | 6 | 12 | 173 | 55 | 95 | 40 | 43 | GREY |
| BT102 | Child | 1 | 11 | 22 | 348 | 97 | 180 | 80 | 85 | GREY |
| BT103 | Small adult | 1 | 16 | 28 | 458 | 130 | 230 | 110 | 120 | GREY |
| BT104 | Adult | 1 | 26 | 33 | 528 | 149 | 260 | 129 | 70 | GREY |
| BT105 | Adult large | 1 | 33 | 41 | 633 | 174 | 330 | 150 | 75 | GREY |
| CT104 | Lower limb cuff adult | 1 | 39 | 55 | 828 | 214 | 440 | 190 | 145 | GREY |
| BT201 | New born | 2 | 6 | 12 | 173 | 55 | 95 | 40 | 43 | GREY |
| BT202 | Child | 2 | 11 | 22 | 348 | 97 | 180 | 80 | 85 | GREY |
| BT203 | Small adult | 2 | 16 | 28 | 458 | 130 | 230 | 110 | 120 | GREY |
| BT204 | Adult | 2 | 26 | 33 | 528 | 149 | 260 | 129 | 70 | GREY |
| BT205 | Adult large | 2 | 33 | 41 | 633 | 174 | 330 | 150 | 75 | GREY |
| CT204 | Lower limb cuff adult | 2 | 39 | 55 | 828 | 214 | 440 | 190 | 145 | GREY |







| | | | |
|---|--|--|------------|
| Dessillons & Dutrillaux Zone Industrielle La Tuque Tél : +33 5 53 48 30 66 technique@ddmedical.fr | 47240 CASTELCUJER Fax : +33 5 53 47 24 64 www.ddmedical.fr | Date of first EC marking certificate : 2004 | Page 2 / 8 |
|---|--|--|------------|

| | | |
|---|--|--|
| A | ↓ Index | Position indication index |
| B |  | Product Brand |
| C | LOT BT020214 | Number of batch |
| D | Child / Enfant REF BT02 | Name and product number |
| E | 11 - 22 CM | Range of the product |
| F | Made in FRANCE €€ | Original product marking |
| G | ← Min ————— Max →  | Validity range of product use |
| H | Artery / Artère Max  | Indication for placement of the artery |
| I |  | Machine wash at 60 ° C |
| J |  | Do not use bleach |
| K |  | Do not iron |
| L | Latex | Latex free |
| M |  | Can be dry cleaned using any solvent |
| N | / 1 Tube | Indication of the number of tube |







A label indicates whether the device has one or two tubes (/ or tube 1/2 tubes).









Connectors for non invasive blood pressure cuffs

| Screw connection | | | |
|-----------------------------|---|-----------------------------|---|
| Reference | Design | Material | Designation |
| A11306 |  | Acetal | Male connector screw connection |
| A11307 |  | Chrome-plated brass | Female connector screw connection |
| A11350 |  | Acetal/ Chrome-plated brass | Male/Female screw connection |
| Quarter-turn connection | | | |
| Reference | Design | Material | Designation |
| A11301 |  | Acetal | Female connector quarter-turn connection |
| A11300 |  | Acetal | Female connector quarter-turn connection with sealing valve |
| A11308 |  | Acetal /EPDM | Male connector quarter-turn connection with O-ring |
| A11351 |  | Acetal/EPDM | Male / female with quarter-turn connection |
| Connection with clip system | | | |
| Reference | Design | Material | Designation |
| A11305 |  | Chrome-plated brass | Female connector with clip system |
| A11313 |  | Chrome-plated brass | Male connector with clip system |
| A11352 |  | Chrome-plated brass | Male / female with clip system |
| A11304 |  | Nylon | Male connector with clip system |
| A11353 |  | Chrome-plated brass /Nylon | Male / female with clip system |


Standard Luer lock connection

| Reference | Design | Material | Designation |
|-----------|--|---------------------|---|
| A11312 |  | Chrome-plated brass | Female connector with luer lock connection |
| A11311 |  | Chrome-plated brass | Male connector with luer lock male connection |
| A11354 |  | Chrome-plated brass | Male / female connector luer lock |
| A11303 |  | Nylon | Female connector luer lock female connection |
| A11302 |  | Nylon | Male connector luer lock male connection |
| A11355 |  | Nylon/ Nylon | Male / female connector luer lock |

Luer lock connection




| Reference | Design | Material | Designation |
|-----------|---|--|--|
| A11050 |  | Nylon | Male luer lock |
| A11356 |  | Nylon/ Nylon | Male / Female luer lock connection |
| A11357 |  | Chrome-plated brass /Nylon | Male / Female connector luer lock /Nylon |
| A11400 |  | Chrome-plated brass | Male luer lock male connection |
| A11358 |  | Chrome-plated brass/ Chrome-plated brass | Male / Female connector luer lock |
| A11359 |  | Nylon/ Chrome-plated brass | Male / Female connector luer lock |

Luer lock Colin connection

| Reference | Design | Material | Designation |
|-----------|---|---------------------|--------------------|
| A11051 |  | Chrome-plated brass | Luer femelle Colin |

Protective caps for Luer lock connection

| Reference | Design | Material | Designation |
|-----------|--------|----------|-------------|
|-----------|--------|----------|-------------|

| | | | |
|--------|---|-----|--|
| A45000 |  | PVC | Protective caps for female luer lock connection |
| A11360 |  | PVC | Female connector luer lock female connection + stopper cap for cleaning and disinfection |
| A11631 |  | PVC | Female connector luer lock female connection + stopper cap for cleaning and disinfection |

3 - Use of the device – Operating instructions

- Select a cuff adapted to the morphology of the patient. The size and range of use are indicated on the cuff.
- Position the center of the bladder pocket represented by a heart on the axis of the brachial artery.
- Roll the cuff around the limb.
- Once the cuff positioned and ready to use, the INDEX must be between the Min and Max (axis of the heart). If INDEX exceeds the Max index located at the axis of the heart, replace the cuff by the model of the larger size. If INDEX exceeds Min index located at the tip of the range, replace the cuff by the model of smaller size.
- When the cuff and INDEX are positioned correctly the procedure for measuring blood pressure can begin.
- Correctly roll the cuff around the limb and ensure its continuation with the fastening system provided.
- Place the tips of the lyre of the stethoscope in your ears, the microphone in the listening position, just below the cuff, on the passage of the artery.
- Ensure that the vent screw of the pear is in the closed position 30 and to inflate or 40 mmHg above the estimated pressure. Blood no longer circulates at the microphone, no sound is heard with a stethoscope.
- Loosen the vent screw the bulb to slowly reduce the inflation pressure and allow blood to flow again, the first beats are perceptible stethoscope.
- Read this pressure gauge reading now, this value corresponds to the maximum pressure or systolic pressure.
- Continue decompression beats levied stethoscope amplify then fade to become inaudible.
- Read this moment the pressure indicated, this reading corresponds to the minimum pressure or diastolic pressure.
- Unscrew the knob further relief to completely purge the remaining air in the bladder.
- Remove the cuff.

DDM

DESSILLONS & DUTRILLAUX - Since 1949

TECHNICAL CHARACTERISTIC


CT 036

Non-invasive blood pressure cuff
made in "Soft Cuff®" fabric

Index 7

4 - Identification and characteristics of different subsets

| Identifier | Designation | Characteristics |
|------------|------------------|--|
| A | Fabric enveloppe | Coated fabric + antibacterial treatment |
| B | Bladder | PU coated fabric - High frequency welded |
| C | Velcro closure | Loop sewn |
| D | Velcro closure | Hook sewn |



Dessillons & Dutrillaux

Zone Industrielle La Tuque
Tél : +33 5 53 48 30 66
technique@ddmedical.fr

47240 CASTELCULIER
Fax : +33 5 53 47 24 64
www.ddmedical.fr

Date of first EC marking certificate :
2004

Page 7 / 8

5 – Cleaning and disinfection

1) Dry cleaning :

- A) Spray a foam detergent, disinfectant, bactericidal on the arm cuff, taking care to distribute the product well on the surface.
- B) Leave 15 minutes.
- C) Do not rinse.
- D) Product recommended: Surface'Safe / Anios.

2) Cleaning by immersion (after removing the bag) :

- A) Remove the arm cuff bag
- B) Clean the arm cuff with water or with washing machine. Washing machine at 40 ° (60 ° being the maximum), without a spin cycle. The number of washing machine and the washing temperature decreases the life of the cuff.
- C) Immerse the arm cuff assembly in a decontaminating solution.
- D) Immersion time : 15 minutes.
Rinse thoroughly with water, renewing the dipping bath after each use. Product recommended : Aniosyme DD1 manufactured by Anios

6 – Packaging

| Reference | Designation | Packaging unit |
|-----------|-----------------------|----------------|
| BT101 | New born | 1 |
| BT102 | Child | 1 |
| BT103 | Small adult | 1 |
| BT104 | Adult | 1 |
| BT105 | Adult large | 1 |
| CT104 | Lower limb cuff adult | 1 |
| BT201 | New born | 1 |
| BT202 | Child | 1 |
| BT203 | Small adult | 1 |
| BT204 | Adult | 1 |
| BT205 | Adult large | 1 |
| CT204 | Lower limb cuff adult | 1 |

8 – Warranty

This warranty provides assurance for the customer who purchases a D & D product that should the product fail to function to D & D published specifications during the term of this warranty, will either replace or repair.

The guarantee period is 1 (one) year from the date of purchase. The product must be used in accordance with its labeling and may not be altered or subjected to misuse, abuse, accident or improper handling.

7 – Storage

| Type of packaging | Storage area | Temperature | Humidity | Atmospheric pressure |
|--------------------|-----------------|---------------|------------|----------------------|
| Original packaging | Ventilated area | -10° à 40 ° c | 30 to 40 % | 500 to 1060 hpA |