

MADE IN FRANCE 

Regulatory information


CE Marking	Index of classification	PHTALATES LATEX
CE	I	

Cuff are used for blood pressure with sphygmomanometer. Arm cuff are made of polyurethane coated fabric and they are immersible.


Flexible and robust design, it's possible to clean in the washing machine and decontaminate them by immersion with the obligation to close the connection tubes.

Reliability: An Index and a range "Min-Max" can quickly validate if the product is adapted to the size of the patient.












1 – Device identification




*Cuff EasyCuff
(BM104)*



*Cuff EasyCuff with A11313 and cap
(BM104A11313)*

A		Indication of the position of the index
B		Product Brand
C		Name and product number
D		Original product marking
E		Range of product use
F		Date of manufacture
G		See operating instructions
H		Latex free
I		Manufacturing Lot Number
J		Validity range of product use
K		Indication for placement of the artery

2 - References



Reference		Designation	Output	Circumference in cm		Dimensions in mm					Tube color
Without connecteur	With connector and cap			Min.	Max.	A	B	C	D	E	
BM101	BM101A11313	Newborn size cuff	1	6	12	170	40	95	47,5	43	WHITE
BM102	BM102A11313	Child size cuff	1	11	22	338	80	180	90	85	PINK
BM103	BM103A11313	Small adult size cuff	1	16	28	423	110	230	115	120	PARMA
BM104	BM104A11313	Adult size cuff	1	26	33	516	129	256	129	70	BLUE SKY
BM105	BM105A11313	Large adult size cuff	1	33	41	616,8	150	330	165	80	TURQUOISE
CM104	CM104A11313	Adult size lower limb cuff	1	39	55	816,5	190	440	220	160	PURPLE
BM201	BM201A11313	Newborn size cuff	2	6	12	170	40	95	47,5	43	WHITE
BM202	BM202A11313	Child size cuff	2	11	22	338	80	180	90	85	PINK
BM203	BM203A11313	Small adult size cuff	2	16	28	423	110	230	115	120	PARMA
BM204	BM204A11313	Adult size cuff	2	26	33	516	129	256	129	70	BLUE SKY
BM205	BM205A11313	Large adult size cuff	2	33	41	616,8	150	330	165	80	TURQUOISE
CM204	CM204A11313	Adult size lower limb cuff	2	39	55	816,5	190	440	220	160	PURPLE

3 – Operating instructions

- Select a cuff adapted to the morphology of the patient. The size and range of use are indicated on the cuff.
- Position the center of the bladder pocket represented by a heart on the axis of the brachial artery.
- Roll the cuff around the limb.
- Once the cuff positioned and ready to use, the INDEX must be between the Min and Max (axis of the heart). If INDEX exceeds the Max index located at the axis of the heart, replace the cuff by the model of the larger size. If INDEX exceeds Min index located at the tip of the range, replace the cuff by the model of smaller size.
- When the cuff and INDEX are positioned correctly the procedure for measuring blood pressure can begin.
- Correctly roll the cuff around the limb and ensure its continuation with the fastening system provided.
- Place the tips of the lyre of the stethoscope in your ears, the microphone in the listening position, just below the cuff, on the passage of the artery.
- Ensure that the vent screw of the pear is in the closed position 30 and to inflate or 40 mmHg above the estimated pressure. Blood no longer circulates at the microphone, no sound is heard with a stethoscope.
- Loosen the vent screw the bulb to slowly reduce the inflation pressure and allow blood to flow again, the first beats are perceptible stethoscope.
- Read this pressure gauge reading now, this value corresponds to the maximum pressure or systolic pressure.
- Continue decompression beats levied stethoscope amplify then fade to become inaudible.
- Read this moment the pressure indicated, this reading corresponds to the minimum pressure or diastolic pressure.
- Unscrew the knob further relief to completely purge the remaining air in the bladder.
- Remove the cuff.

4 - Maintenance:

- Under no circumstances should liquid enter the pneumatic part of the cuff or shorts, as this could cause partial or total damage to the electro-medical device or the hand-held measuring device to which the cuff or shorts are connected.

In order to avoid any risk, we recommend the D&D sealing systems:

- Plugs, Clamps, D&D Quick Valve System.
- Under no circumstances should a tube be clamped with pliers or other unsuitable instruments.
- Under no circumstances should the cuff or the shorts be manipulated with any sharp object.
- Under no circumstances should the cuff or shorts be inflated without being rolled up and closed with the two Velcro parts on a limb or round tube without edges and sharp parts.

5 - Cleaning:

EasyCuff Cuffs and Shorts are not sterilizable.

Reference	Dry cleaning	Cleaning by immersion
BM101	X	
BM102	X	
BM103	X	
BM104	X	
BM105	X	
CM104	X	
BM201	X	
BM202	X	
BM203	X	
BM204	X	
BM205	X	
CM204	X	
BM101A11313	X	X
BM102A11313	X	X
BM103A11313	X	X
BM104A11313	X	X
BM105A11313	X	X
CM104A11313	X	X
BM201A11313	X	X
BM202A11313	X	X
BM203A11313	X	X
BM204A11313	X	X
BM205A11313	X	X
CM204A11313	X	X

1) Dry cleaning:

- A) Spray a detergent, disinfectant, bactericidal foam on the arm cuff, taking care to distribute the product well on the surface.
- B) Leave for 15 minutes.
- C) Do not rinse.
- D) Product recommended: Surface'Safe / Anios.

2) Cleaning by immersion :

- A) Seal the tube or tubes by means of the closure system provided:
- Connector with cap: A11313
- B) Immerse the tube and arm cuff assembly in a decontaminating solution.
- C) Immersion time : 15 minutes.
- D) Rinse thoroughly with water, renewing the dipping bath after each use. Product recommended: Aniosyme DD1 manufactured by Anios.

6 - Storage

Type of packaging	Storage area	Temperature	Humidity	Atmospheric pressure
Original packaging	Ventilated area	-10° to 40 ° c	30 to 40 %	500 to 1060 hpA

 <p>DDM DESSILLONS & DUTRILLAUX</p>	INSTRUCTIONS FOR USE	INS 104
	Immersible arm-cuff "Easy Cuff"	Index 7

7 - Warranty

This warranty provides assurance for the customer who purchases a D & D product that should the product fail to function to D & D published specifications during the term of this warranty, will either replace or repair.

The guarantee period is 1 (one) year from the date of purchase.

The product must be used in accordance with its labeling and may not be altered or subjected to misuse, abuse, accident or improper handling.