

	<b>TECHNICAL CHARACTERISTIC</b>	CT 009
	Immersible arm-cuff "Easy Cuff"	Index 10

**MADE IN FRANCE** 

**Regulatory information**


CE Marking	Index of classification	<b>PHTALATES</b>	<b>LATEX</b>
CE	I		

Arm cuff are polyurethane coated fabric and are immersible.


Flexible and robust design, it's possible to clean in the washing machine and decontaminate them by immersion with the obligation to close the connection tubes.

Reliability: An Index and a range "Min-Max" can quickly validate if the product is adapted to the size of the patient.

**1 - Device identification**














*Cuff EasyCuff  
(BM104)*




*Cuff EasyCuff with A11313 and cap  
(BM104A11313)*

<b>Dessillons &amp; Dutrillaux</b>		Date of first EC marking certificate : <b>2012</b>	Page 1 / 10
Zone Industrielle La Tuque Tél : +33 5 53 48 30 66 technique@ddmedical.fr	47240 CASTELCULIER Fax : +33 5 53 47 24 64 www.ddmedical.fr		

A		Indication of the position of the index
B		Product Brand
C		Name and product number
D		Original product marking
E		Range of product use
F		Date of manufacture
G		See operating instructions
H		Latex free
I		Manufacturing Lot Number
J		Validity range of product use
K		Indication for placement of the artery

## 2 – References

According to EN 1060-1 EC, the optimal size of the inflation bag must be at least equal in width to 40% of the circumference of the arm, and a length between 80 and 100% of the circumference of the arm (see markings shown on cuffs Soft Cuff).



Reference		Designation	Output	Circumference in cm		Dimensions in mm					Tube color
Without connecteur	With connector and cap			Min.	Max.	A	B	C	D	E	
BM101	BM101A11313	Newborn size cuff	1	6	12	170	40	95	47,5	43	WHITE
BM102	BM102A11313	Child size cuff	1	11	22	338	80	180	90	85	PINK
BM103	BM103A11313	Small adult size cuff	1	16	28	423	110	230	115	120	PARMA
BM104	BM104A11313	Adult size cuff	1	26	33	516	129	256	129	70	BLUE SKY
BM105	BM105A11313	Large adult size cuff	1	33	41	616,8	150	330	165	80	TURQUOISE
CM104	CM104A11313	Adult size lower limb cuff	1	39	55	816,5	190	440	220	160	PURPLE
BM201	BM201A11313	Newborn size cuff	2	6	12	170	40	95	47,5	43	WHITE
BM202	BM202A11313	Child size cuff	2	11	22	338	80	180	90	85	PINK
BM203	BM203A11313	Small adult size cuff	2	16	28	423	110	230	115	120	PARMA
BM204	BM204A11313	Adult size cuff	2	26	33	516	129	256	129	70	BLUE SKY
BM205	BM205A11313	Large adult size cuff	2	33	41	616,8	150	330	165	80	TURQUOISE
CM204	CM204A11313	Adult size lower limb cuff	2	39	55	816,5	190	440	220	160	PURPLE

 <p><b>DDM</b> DESSILLONS &amp; DUTRILLAUX - Since 1949</p>	<b>TECHNICAL CHARACTERISTIC</b>	CT 009
	Immersible arm-cuff "Easy Cuff"	Index 10

### 3 - Use of the device – Operating instructions

- Select a cuff adapted to the morphology of the patient. The size and range of use are indicated on the cuff.
- Position the center of the bladder pocket represented by a heart on the axis of the brachial artery.
- Roll the cuff around the limb.
- Once the cuff positioned and ready to use, the INDEX must be between the Min and Max (axis of the heart). If INDEX exceeds the Max index located at the axis of the heart, replace the cuff by the model of the larger size. If INDEX exceeds Min index located at the tip of the range, replace the cuff by the model of smaller size.
- When the cuff and INDEX are positioned correctly the procedure for measuring blood pressure can begin.
- Correctly roll the cuff around the limb and ensure its continuation with the fastening system provided.
- Place the tips of the lyre of the stethoscope in your ears, the microphone in the listening position, just below the cuff, on the passage of the artery.
- Ensure that the vent screw of the pear is in the closed position 30 and to inflate or 40 mmHg above the estimated pressure. Blood no longer circulates at the microphone; no sound is heard with a stethoscope.
- Loosen the vent screw the bulb to slowly reduce the inflation pressure and allow blood to flow again, the first beats are perceptible stethoscope.
- Read this pressure gauge reading now, this value corresponds to the maximum pressure or systolic pressure.
- Continue decompression beats levied stethoscope amplify then fade to become inaudible.
- Read this moment the pressure indicated, this reading corresponds to the minimum pressure or diastolic pressure.
- Unscrew the knob further relief to completely purge the remaining air in the bladder.
- Remove the cuff.

<b>Dessillons &amp; Dutrillaux</b>		Date of first EC marking certificate : <b>2012</b>	Page 4 / 10
Zone Industrielle La Tuque Tél : +33 5 53 48 30 66 technique@ddmedical.fr	47240 CASTELCULIER Fax : +33 5 53 47 24 64 www.ddmedical.fr		

**DDM**

DESSILLONS &amp; DUTRILLAUX - Since 1949

**TECHNICAL CHARACTERISTIC**

CT 009

Immersible arm-cuff "Easy Cuff"

Index 10

**4 - Identification and characteristics of different subsets**

Identifier	Designation	Characteristics
A	Immersible monobloc assembly cuff adult	PU coated fabric - high-frequency welding
B	Velcro	Loop Astrakan high-frequency welding
C	Angled connecting tube	Molded PU, welding HF, angle 90 °
D	Connecting hose	Compound vinylique phtalates free 58 shore A
E	Female connector luer lock female connection + stopper cap for cleaning and disinfection	PVC / Nylon
F	Male fitting with clip connection (reference A11313) with cap	Chrome plated brass

Dessillons &amp; Dutrillaux

Zone Industrielle La Tuque

Tél : +33 5 53 48 30 66

technique@ddmedical.fr


47240 CASTELCULIER

Fax : +33 5 53 47 24 64


www.ddmedical.fr

Date of first EC marking certificate :  
2012

Page 5 / 10



Cuff EasyCuff



Cuff EasyCuff with A11313 andcap

## 5 – Cleaning and disinfection

*EasyCuff Cuffs and Shorts are not sterilizable.*

Reference	Dry cleaning	Cleaning by immersion
BM101	X	
BM102	X	
BM103	X	
BM104	X	
BM105	X	
CM104	X	
BM201	X	
BM202	X	
BM203	X	
BM204	X	
BM205	X	
CM204	X	
BM101A11313	X	X
BM102A11313	X	X
BM103A11313	X	X
BM104A11313	X	X
BM105A11313	X	X
CM104A11313	X	X
BM201A11313	X	X
BM202A11313	X	X
BM203A11313	X	X
BM204A11313	X	X
BM205A11313	X	X
CM204A11313	X	X

	<b>TECHNICAL CHARACTERISTIC</b>	CT 009
	Immersible arm-cuff "Easy Cuff"	Index 10

**1) Dry cleaning:**

- A) Spray a detergent, disinfectant, bactericidal foam on the arm cuff, taking care to distribute the product well on the surface.
- B) Leave for 15 minutes.
- C) Do not rinse.
- D) Product recommended: Surface'Safe / Anios.

**2) Cleaning by immersion :**

- A) Seal the tube or tubes by means of the closure system provided:
  - Connector with cap: A11313
- B) Immerse the tube and arm cuff assembly in a decontaminating solution.
- C) Immersion time : 15 minutes.
- D) Rinse thoroughly with water, renewing the dipping bath after each use. Product recommended: Aniosyme DD1 manufactured by Anios.

<b>Dessillons &amp; Dutrillaux</b>		Date of first EC marking certificate : <b>2012</b>	<b>Page 8 / 10</b>
<b>Zone Industrielle La Tuque</b> Tél : +33 5 53 48 30 66 <b>technique@ddmedical.fr</b>	<b>47240 CASTELCULIER</b> Fax : +33 5 53 47 24 64 <b>www.ddmedical.fr</b>		



### 6 – Packaging

Reference		Designation	Packaging unit
Without connector	With connector and cap		
BM101	BM101A11313	Newborn size cuff	1
BM102	BM102A11313	Child size cuff	1
BM103	BM103A11313	Small adult size cuff	1
BM104	BM104A11313	Adult size cuff	1
BM105	BM105A11313	Large adult size cuff	1
CM104	CM104A11313	Adult size lower limb	1
BM201	BM201A11313	Newborn size cuff	1
BM202	BM202A11313	Child size cuff	1
BM203	BM203A11313	Small adult size cuff	1
BM204	BM204A11313	Adult size cuff	1
BM205	BM205A11313	Large adult size cuff	1
CM204	CM204A11313	Adult size lower limb	1

### 7 – Storage

Type of packaging	Storage area	Temperature	Humidity	Atmospheric pressure
Original packaging	Ventilated area	-10° to 40 ° c	30 to 40 %	500 to 1060 hpA

 <p>DESSILLONS &amp; DUTRILLAUX - Since 1949</p>	<p><b>TECHNICAL CHARACTERISTIC</b></p>	<p>CT 009</p>
	<p>Immersible arm-cuff "Easy Cuff"</p>	<p>Index 10</p>

## 8 – Warranty

This warranty provides assurance for the customer who purchases a D & D product that should the product fail to function to D & D published specifications during the term of this warranty, will either replace or repair.

The guarantee period is 1 (one) year from the date of purchase.

The product must be used in accordance with its labeling and may not be altered or subjected to misuse, abuse, accident or improper handling.

<p>Dessillons &amp; Dutrillaux</p>		<p>Date of first EC marking certificate : 2012</p>	<p>Page 10 / 10</p>
<p>Zone Industrielle La Tuque Tél : +33 5 53 48 30 66 technique@ddmedical.fr</p>	<p>47240 CASTELCULIER Fax : +33 5 53 47 24 64 www.ddmedical.fr</p>		