

MADE IN FRANCE 

Regulatory information

| CE Marking | Index of classification | Notified body | PHTALATES LATEX |
|------------|-------------------------|---------------|------------------------|
| CE 0459 | Ila | G-MED | |



Technical characteristics Easy Pump G10903

| | |
|----------------------|--|
| STORAGE CONDITIONS | T : -5° à 40°C, Humidity : 20- 80 % |
| CONDITIONS OF USE | T : 5° à 40°C, Humidity : 20 – 80 % Maximum altitude 2000 m |
| CASE MATERIALS | ABS (Acrylonitrile Butadiene Styrene) |
| Density | 1.112 g/ cm ³ |
| Shock resistance | > 25 KJ / m ² |
| UL flammability test | UL94 V-0 |
| Heat stability | 85°C |
| Hardness | 78 |
| DIMENSIONS | CASE |
| Height (mm) | 250 |
| Width (mm) | 300 |
| Depth (mm) | 190 |
| | SCREEN |
| Width (mm) | 120 |
| Height (mm) | 90 |
| WEIGHT kg | 4.2 KGS |
| LINE VOLTAGE | 100 - 240 V AC |
| LINE FREQUENCY | 50 – 60 Hz |
| POWER PLUG | 60 VA |
| FUSE | |
| Input power card | F2AH/250V |
| Output power card | F1,6AH/15V |
| Battery | F8AH/12V |
| PULSED POWER SUPPLY | Entry : 90-264 Vac, 47-63 Hz, 1.8 – 1 A Exit : 15 Vdc, 4A |
| BATTERY | Rechargeable |
| Type | Pb |
| Voltage | 12 Volts |
| Ampere | 4000 mAh |
| Charging time | 8 h |
| Autonomy | 10 h |
| PUMP | |
| Type | Membrane pump |
| Pump flow | 4.6 l/mn |

DDM

ISO 13485



DESSILLONS & DUTRILLAUX

TECHNICAL CHARACTERISTICS

CT 041

Pneumatic tourniquet G10903 Easy Pump
dual 2 tubes

Index 4

| | | |
|--|-----------|---------------------------------------|
| PRESSION | | |
| Type | | 1.7 bar |
| Service | | 0 to 600 mmHg |
| Setting range | | ± 1 mm Hg |
| Setting precision | | ± 5 mm Hg |
| | | An audible and visual alarm |
| | | 2 |
| TIMER | | |
| Units | | minutes |
| Alarm | | Programmable audible and visual alarm |
| CONNECTION | | |
| | Pneumatic | 2 female coupler CPC type |
| | Electric | Connector CEE22 |
| MAXIMUM SOUND LEVEL WHILE FUNCTIONING | | 52 dB |

Dessillons & Dutrillaux

Zone Industrielle La Tuque

47240 CASTELCUJER

Tél : +33 5 53 48 30 66

Fax : +33 5 53 47 24 64

technique@ddmedical.fr

www.ddmedical.fr

Date of first EC marking certificate :
1998

Page 3 / 3