

**DDM**

ISO 13485



DESSILLONS &amp; DUTRILLAUX

**TECHNICAL CHARACTERISTICS**

CT 016

Stethoscope Bi-Sonal  
with dual-head chestpiece

Index 6

**Regulatory information**

CE Marking	Index of classification	<del><b>PHTALATES LATEX</b></del>
CE	I	

Stethoscope with dual-head chestpiece made of stainless steel equipped with a non-chill ring and flexible ear-tips.

Binaural is made of chrome plated brass.

**1 – Device identification****2 – References**

Reference	Designation
S20019	Stethoscope dual-head chestpiece Ø 17 mm Newborn
S20011	Stethoscope dual-head chestpiece Ø 34 mm Pediatric
S20010	Stethoscope dual-head chestpiece Ø 45 mm Adult

Dessillons &amp; Dutrillaux

Zone Industrielle La Tuque

Tél : +33 5 53 48 30 66

technique@ddmedical.fr

47240 CASTELCULIER

Fax : +33 5 53 47 24 64

www.ddmedical.fr

Date of first EC marking certificate :  
2004


Page 1 / 3

### 3 – Use of the device - Operating instructions

- ⇒ Heart auscultation                   ■ ■ □
- ⇒ Vascular auscultation               ■ ■ □
- ⇒ Pulmonary auscultation           ■ ■ □
- ⇒ Blood pressure monitoring       ■ ■ □

### 4 - Identification and characteristics of different subsets :

Identifiant	Reference	Designation	Characteristic
A	A10360	Dual-head chestpiece Ø 17 mm	Stainless steel
A	A10109	Dual-head chestpiece Ø 34 mm	Stainless steel
A	A10108	Dual-head chestpiece Ø 45 mm	Stainless steel
B	A10361	Diaphragm Ø 17 mm	Terphane
B	A10351	Diaphragm Ø 34 mm	Terphane
B	A10106	Diaphragm Ø 45 mm	Terphane
C	A10103	Ear tips	Flexible moulded PVC
D	A10104	Binaural (with ear tips)	Chrome plated-brass tube + PVC
E	A10105	Tubing	Moulded PVC
F	A10350	Non-chill ring	PVC



### 5 – Cleaning and disinfection

Spray disinfectant solutions. Do not use detergent, do not immerse.

### 6 - Packaging

Reference	Packaging unit
S20019	1
S20011	1
S20010	1

### 7 - Storage

Type of packaging	Storage area	Temperature	Humidity	Atmospheric pressure
Original packaging	Ventilated area	-10° to 40 ° c	30 to 40 %	500 to 1060 hpA



## TECHNICAL CHARACTERISTICS

CT 016

Stethoscope Bi-Sonal  
with dual-head chestpiece

Index 6

### 8 - Warranty

This warranty provides assurance for the customer who purchases a D & D product that should the product fail to function to D & D published specifications during the term of this warranty, will either replace or repair.

The duration of this warranty is 12 months.

The product must be used in accordance with its labeling and may not be altered or subjected to misuse, abuse, accident or improper handling.

Dessillons & Dutrillaux

Zone Industrielle La Tuque

Tél : +33 5 53 48 30 66

technique@ddmedical.fr

47240 CASTELCULIER

Fax : +33 5 53 47 24 64

www.ddmedical.fr

Date of first EC marking certificate :  
2004

Page 3 / 3