

### Regulatory information

| CE Marking | Index of classification | <del>LATEX</del> |
|------------|-------------------------|------------------|
| CE         | I                       |                  |

Stethoscope with single head chestpiece made of stainless steel equipped with a non-chill ring and flexible ear-tips.

Binaural is made of chrome plated brass.

### 1 – Device identification



### 2 – References


| Reference | Designation  |
|-----------|--|
| S20032    | Stethoscope single head chestpiece Ø 34 mm Pediatric |
| S20008    | Stethoscope single head chestpiece Ø 45 mm Adult     |

### 3 – Use of the device - Operating instructions

- ⇒ Heart auscultation                   ■ ■ □
- ⇒ Vascular auscultation               ■ ■ □
- ⇒ Pulmonary auscultation           ■ ■ □
- ⇒ Blood pressure monitoring       ■ ■ □

### 4 - Identification and characteristics of different subsets :

| Identifier | Reference | Designation                    | Characteristic                      |
|------------|-----------|--------------------------------|-------------------------------------|
| A          | A10119    | Single head chestpiece Ø 34 mm | Stainless steel                     |
| A          | A10102    | Single head chestpiece Ø 45 mm | Stainless steel                     |
| B          | A10124    | Diaphragm Ø 34 mm              | Terphane                            |
| B          | A10106    | Diaphragm Ø 45 mm              | Terphane                            |
| C          | A10103    | Ear tips                       | Flexible moulded PVC                |
| D          | A10104    | Binaural (with ear tips)       | Chrome plated-brass tube tube + PVC |
| E          | A10105    | Tubing                         | Moulded PVC                         |
| F          | A10350    | Non-chill ring                 | PVC                                 |



### 5 – Cleaning and disinfection

Spray disinfectant solutions. Do not use detergent, do not immerse.

### 6 - Packaging

| Reference | Packaging unit |
|-----------|----------------|
| S20032    | 1              |
| S20008    | 1              |

### 7 - Storage

| Type of packaging  | Storage area    | Temperature    | Humidity   | Atmospheric pressure |
|--------------------|-----------------|----------------|------------|----------------------|
| Original packaging | Ventilated area | -10° to 40 ° c | 30 to 40 % | 500 to 1060 hpA      |



## TECHNICAL CHARACTERISTICS

CT 017

Stethoscope SONAL with single head  
chestpiece

Index 6

### 8 - Warranty

This warranty provides assurance for the customer who purchases a D & D product that should the product fail to function to D & D published specifications during the term of this warranty, will either replace or repair.

The duration of this warranty is 12 months.

The product must be used in accordance with its labeling and may not be altered or subjected to misuse, abuse, accident or improper handling.

Dessillons & Dutrillaux

Zone Industrielle La Tuque

Tél : +33 5 53 48 30 66

technique@ddmedical.fr

47240 CASTELCULIER

Fax : +33 5 53 47 24 64

www.ddmedical.fr

Date of first EC marking certificate :  
2004

Page 3 / 3