

### Regulatory information

CE Marking	Index of classification	<del>LATEX</del>
CE	I	

Stethoscope with single head chestpiece made of aluminium equipped with a non-chill ring and flexible ear-tips.  
 Binaural is made of chrome plated brass.

### 1 – Device identification



### 2 – References

Reference	Designation
S20034	Stethoscope single head chestpiece Ø 34 mm Pediatric
S20003	Stethoscope single head chestpiece Ø 45 mm Adult

Dessillons & Dutrillaux

Zone Industrielle La Tuque  
 Tél : +33 5 53 48 30 66  
 technique@ddmedical.fr

47240 CASTELCULIER  
 Fax : +33 5 53 47 24 64  
 www.ddmedical.fr

Date of first EC marking certificate :  
 2004


Page 1 / 3

### 3 – Use of the device - Operating instructions

- ⇒ Heart auscultation                   ■ □ □
- ⇒ Vascular auscultation               ■ □ □
- ⇒ Pulmonary auscultation           ■ □ □
- ⇒ Blood pressure monitoring       ■ □ □

### 4 - Identification and characteristics of different subsets :

Identifier	Reference	Designation	Characteristic
A	A10122	Single head chestpiece Ø 34 mm	Aluminium
A	A12107	Single head chestpiece Ø 45 mm	Aluminium
B	A10123	Diaphragm Ø 34 mm	Terphane
B	A10106	Diaphragm Ø 45 mm	Terphane
C	A10103	Ear tips	Flexible moulded PVC
D	A12104	Binaural (with ear tips)	Chrome plated-brass tube + PVC
E	A12105	Tubing	Moulded PVC
F	A12401	Non-chill ring	PVC



### 5 – Cleaning and disinfection

Spray disinfectant solutions. Do not use detergent, do not immerse.

### 6 - Packaging

Reference	Packaging unit
S20034	1
S20003	1

### 7 - Storage

Type of packaging	Storage area	Temperature	Humidity	Atmospheric pressure
Original packaging	Ventilated area	-10° to 40 ° c	30 to 40 %	500 to 1060 hpA

 <p><b>DDM</b> DESSILLONS &amp; DUTRILLAUX</p>	<b>TECHNICAL CHARACTERISTICS</b>	CT 018
	Stethoscope SOLO with single head chestpiece	Index 6

## 8 - Warranty

This warranty provides assurance for the customer who purchases a D & D product that should the product fail to function to D & D published specifications during the term of this warranty, will either replace or repair.

The duration of this warranty is 12 months.

The product must be used in accordance with its labeling and may not be altered or subjected to misuse, abuse, accident or improper handling.

<p align="center"><b>Dessillons &amp; Dutrillaux</b></p>		<p>Date of first EC marking certificate : 2004</p>	<p align="right">Page 3 / 3</p>
<p>Zone Industrielle La Tuque Tél : +33 5 53 48 30 66 technique@ddmedical.fr</p>	<p>47240 CASTELCULIER Fax : +33 5 53 47 24 64 www.ddmedical.fr</p>		