

Regulatory information

CE Marking	Index of classification	LATEX
CE	I	

Stethoscope with single head chestpiece made of brushed brass equipped with a non-chill ring and flexible ear-tips.

Binaural is made of chrome plated brass.

1 – Device identification



2 – References


Reference	Designation
S20000	Stethoscope single head chestpiece Ø 45 mm Adult

3 – Use of the device - Operating instructions

- ⇒ Heart auscultation ■ ■ □
- ⇒ Vascular auscultation ■ ■ □
- ⇒ Pulmonary auscultation ■ ■ □
- ⇒ Blood pressure monitoring ■ ■ □

4 - Identification and characteristics of different subsets :

Identifier	Reference	Designation	Characteristic
A	A12407	Single head chestpiece Ø 45 mm	Brushed brass
B	A12406	Diaphragm Ø 45 mm	Terphane
C	A10103	Ear tips	Flexible moulded PVC
D	A10104	Binaural (with ear tips)	Chrome plated-brass tube + PVC
E	A12405	Tubing	Moulded PVC
F	A12401	Non-chill ring	PVC



5 – Cleaning and disinfection

Spray disinfectant solutions. Do not use detergent, do not immerse.

6 - Packaging

Reference	Packaging unit
S20000	1

7 - Storage

Type of packaging	Storage area	Temperature	Humidity	Atmospheric pressure
Original packaging	Ventilated area	-10° to 40 ° c	30 to 40 %	500 to 1060 hpA

8 - Warranty

This warranty provides assurance for the customer who purchases a D & D product that should the product fail to function to D & D published specifications during the term of this warranty, will either replace or repair.

The duration of this warranty is 12 months.

The product must be used in accordance with its labeling and may not be altered or subjected to misuse, abuse, accident or improper handling.