

Regulatory information

CE Marking	Index of classification	LATEX
CE	I	

Stethoscope with dual-head chestpiece made of aluminium with a non-chill ring and flexible ear-tips. Binaural is made of chrome plated brass.

1 – Device identification



2 – References


Reference	Designation
S20033	Stethoscope dual-head chestpiece Ø 17 mm Newborn
S20035	Stethoscope dual-head chestpiece Ø 34 mm Pediatric
S20015	Stethoscope dual-head chestpiece Ø 45 mm Adult

3 – Use of the device - Operating instructions

- ⇒ Heart auscultation ■ □ □
- ⇒ Vascular auscultation ■ □ □
- ⇒ Pulmonary auscultation ■ □ □
- ⇒ Blood pressure monitoring ■ □ □

4 - Identification and characteristics of different subsets :

Identifiant	Reference	Designation	Characteristic
A	A10120	Dual-head chestpiece Ø 17 mm	Aluminium
A	A10121	Dual-head chestpiece Ø 34 mm	Aluminium
A	A12207	Dual-head chestpiece Ø 45 mm	Aluminium
B	A10126	Diaphragm Ø 17 mm	Terphane
B	A10125	Diaphragm Ø 34 mm	Terphane
B	A12206	Diaphragm Ø 45 mm	Terphane
C	A12203	Ear tips	Flexible moulded PVC
D	A12204	Binaural (with ear tips)	Chrome plated-brass tube + PVC
E	A12201	Non-chill ring	PVC



5 – Cleaning and disinfection

Spray disinfectant solutions. Do not use detergent, do not immerse.

6 - Packaging

Reference	Packaging unit
S20033	1
S20035	1
S20015	1

7 - Storage

Type of packaging	Storage area	Temperature	Humidity	Atmospheric pressure
Original packaging	Ventilated area	-10° to 40 ° c	30 to 40 %	500 to 1060 hpA

8 - Warranty

 <p>DDM DESSILLONS & DUTRILLAUX</p>	TECHNICAL CHARACTERISTICS	CT 020
	Stethoscope VEGA with dual-head chestpiece	Index 6

This warranty provides assurance for the customer who purchases a D & D product that should the product fail to function to D & D published specifications during the term of this warranty, will either replace or repair. The duration of this warranty is 12 months. The product must be used in accordance with its labeling and may not be altered or subjected to misuse, abuse, accident or improper handling.

Dessillons & Dutrillaux Zone Industrielle La Tuque Tél : +33 5 53 48 30 66 technique@ddmedical.fr		47240 CASTELCULIER Fax : +33 5 53 47 24 64 www.ddmedical.fr	Date of first EC marking certificate : 2004	Page 3 / 3
---	--	---	--	------------