

### Regulatory information

CE Marking	Index of classification	<del><b>PHTALATES LATEX</b></del>
CE 0120	I	

Stethoscope with dual-head chestpiece made of aluminium and equipped with an aluminium screwed chestpiece and flexible ear-tips.

The binaural is made of chromed-plated brass with dual PVC tubes.

### 1 - Device identification



### 2 – References


Reference	Designation
S20028	Stethoscope dual-head chestpiece Ø 45 mm Adult

### 3 - Use of the device – Operating instructions

- Heart auscultation                   ■ ■ □
- Vascular auscultation               ■ ■ □
- Auscultation of lungs               ■ ■ □
- Blood pressure measurement       ■ ■ □

### 4 - Identification and characteristics of different subsets

Identifier	Reference	Designation	Characteristics
A	A12307	Dual-head chestpiece Ø 45 mm	Aluminium
B	A12306	Diaphragm Ø 45 mm	Terphane
C	A12303	Ear tips	Moulded PVC
D	A12304	Binaural (with ear tips)	Chrome plated brass tube + PVC
E	A12305	Double tubing	Moulded PVC
F	A12301	Chestpiece ring	Aluminium



### 5 – Cleaning and disinfection

Spray disinfectant solutions. Do not use detergent, do not immerse.

### 6 – Packaging

Reference	Packaging unit
S20028	1

### 7 – Storage

Type of packaging	Storage area	Temperature	Humidity	Atmospheric pressure
Original packaging	Ventilated area	-10° to 40 ° c	30 to 40 %	500 to 1060 hpA

## 8 – Warranty

This warranty provides assurance for the customer who purchases a D & D product that should the product fail to function to D & D published specifications during the term of this warranty, will either replace or repair.

The guarantee period is 12 months from the date of purchase. The product must be used in accordance with its labeling and may not be altered or subjected to misuse, abuse, accident or improper handling.